

# Biosensors

**B. Tech.**

**Course No.: EEL 3050**

**L-T-P [C]: 3-0-2 [4]**

***Prof. AJAY AGARWAL***

**ELECTRICAL ENGINEERING**

**IIT JODHPUR**


*Lecture 22 dated 3<sup>rd</sup> Oct. 2024*

## Electrochemical Sensors:

- The **top two** most popular transducers for biosensors are **optical & electrochemical**, followed by piezoelectric & thermal.
- **Electrochemical sensors** are electrochemical transducers, where the **conc. of ions or chemicals** are converted into
  - electrical **voltage** (potentiometric),
  - electrical **current** (amperometric), or
  - electrical **resistance/ conductance** (conductometric).
- If **electrochemical** sensors are **used** together **with** **bioreceptors** (**enzymes** or **antibodies**), they become **electrochemical biosensors**.

## Electrolytic and Electrochemical Cells:


- An **electrolytic cell** decomposes **ionic chemical compounds** by applying **voltage** to its solution.
- a typical electrolytic cell has **two electrodes** (metal rods) **inserted** into a **solution of metal salts** (electrolytes).
- **Electrons are taken** from metal ions at one electrode (**oxidation**) and are **released** to metal ions at the other electrode (**reduction**).
- Whole reaction is called redox (reduction + oxidation) reaction.
- The current that flows between the electrodes depends not only on the voltage that is applied, but also on the electrical properties of the solution.



**An electrolytic cell**


## Electrochemical Cells:

- An **electrochemical cell** is a device that **generates** electrical voltage and current.
- Two **electrodes** are **inserted** into **electrolytes** that are separated by a **salt bridge** (Galvanic cell), or a **semi-permeable membrane** that replaces the function of a salt bridge (Daniell cell)
- The **cells are identical** to each other
- **Each side** of the cells is referred to as a **half-cell**.
- The cells are described as:



Electrochemical cells: galvanic (left) and Daniell (right) cells

The following chemical formulas describe the equilibrium of both metal ions:



Combining these two formulas yields:



Here, forward reaction is spontaneous i.e. the change in Gibbs free energy is negative or  $\Delta G < 0$ .  $\Delta G$  is defined as

$$\Delta G = -n.F.\Delta E \quad \dots \quad (3)$$

where,

n number of electrons

F Faraday constant =  $96,487 \text{ C mol}^{-1}$

$\Delta E$  electrical potential difference (V)

# Biosensors

**B. Tech.**

**Course No.: EEL 3050**

**L-T-P [C]: 3-0-2 [4]**

***Prof. AJAY AGARWAL***


**ELECTRICAL ENGINEERING**

**IIT JODHPUR**

*Lecture 23 dated 16<sup>th</sup> Oct. 2024*


## Electrochemical Cells:

- An **electrochemical cell** is a device that **generates** electrical voltage and current.
- Two **electrodes** are **inserted** into **electrolytes** that are separated by a **salt bridge** (Galvanic cell), or a **semi-permeable membrane** that replaces the function of a salt bridge (Daniell cell)
- The **cells are identical** to each other
- **Each side** of the cells is referred to as a **half-cell**.
- The cells are described as:



Electrochemical cells: galvanic (left) and Daniell (right) cells

The following chemical formulas describe the equilibrium of both metal ions:



Combining these two formulas yields:



Here, forward reaction is spontaneous i.e. the change in Gibbs free energy is negative or  $\Delta G < 0$ .  $\Delta G$  is defined as

$$\Delta G = -n.F.\Delta E \quad \dots \quad (3)$$

where,

n number of electrons

F Faraday constant =  $96,487 \text{ C mol}^{-1}$

$\Delta E$  electrical potential difference (V)



We also know that

$$\Delta G = RT \ln K \quad \dots \quad (4)$$

where,

R gas constant = 8.3145 J K<sup>-1</sup> mol<sup>-1</sup>

T temperature (K)

K equilibrium constant

Combining above Eqs. 3 & 4 gives

$$\Delta E = -RT / nF \times \ln K \quad \dots \quad (5)$$

Applying Eq. 5 to the left-side half-cell of the above electrochemical cell:  $\text{Cu}^{2+} + 2\text{e}^- = \text{Cu}$ .

For this reaction,  $n$  is 2, and  $K$  is defined by the **activity ratio**

$$K = a_{\text{Cu}} / a_{\text{Cu}^{2+}} \quad \dots \quad (6)$$

The activity of a pure liquid or pure solid is 1 ( $a_{\text{Cu}} = 1$ ), and the activity of metal ions can be approximated by their molar concentration ( $a_{\text{Cu}^{2+}} = [\text{Cu}^{2+}]$ ).

$\Delta E$  can also be replaced with  $E - E^\circ$ , where  $E^\circ$  refers to the **standard electrode potential** (reference or ground potential).

Plugging all of these into Eq. 10.5 yields

$$E = E^\circ - \frac{RT}{2F} \times \ln 1 / [\text{Cu}^{2+}]$$

or

$$E = E^\circ + \frac{RT}{2F} \times \ln [\text{Cu}^{2+}] \quad \dots \quad (7)$$

Or more generally,

$$E = E^\circ + \frac{RT}{nF} \times \ln [M^+] \quad \dots \quad (8)$$

which is called the **Nernst equation**. Notice that the **potential (voltage)** is **linearly proportional** to the **logarithm of ionic concentration**.

We can also convert the **natural log** into a **common log** using  $\ln X = 2.303 \log X$ . By plugging in universal constants for  $R$ ,  $F$ , and a room temperature of  $25\text{ }^{\circ}\text{C}$  ( $=273.15\text{ K}$ ) into Eq. 8

$$E = E^{\circ} + 0.059/n \times \log [M^+] \quad \dots \quad (9)$$

For the entire cell,



$$K = a_{\text{Cu}} a_{\text{Zn}^{2+}} / a_{\text{Cu}^{2+}} a_{\text{Zn}} = a_{\text{Zn}^{2+}} / a_{\text{Cu}^{2+}} \approx [\text{Zn}^{2+}] / [\text{Cu}^{2+}] \quad \dots \quad (10)$$

as the activities of pure solids are 1 ( $a_{\text{Cu}} = a_{\text{Zn}} = 1$ ).

Plugging Eq. 10 into Eq. 5 and repeating the above procedure yields

$$E = E^\circ + 0.059/n \times \log [M^+] \quad \dots \quad (9)$$

For the entire cell,



$$K = a_{\text{Cu}} a_{\text{Zn}^{2+}} / a_{\text{Cu}^{2+}} a_{\text{Zn}} = a_{\text{Zn}^{2+}} / a_{\text{Cu}^{2+}} \approx [\text{Zn}^{2+}] / [\text{Cu}^{2+}] \quad \dots \quad (10)$$

as the activities of pure solids are 1 ( $a_{\text{Cu}} = a_{\text{Zn}} = 1$ ). Plugging Eq. 10 into Eq. 5 & repeating the above procedure yields

$$E = E^\circ + 0.059/n \times \log [\text{Zn}^{2+}] / [\text{Cu}^{2+}] \quad \dots \quad (11)$$

For a Galvanic or Daniell cell, the standard electrode potential for the entire reaction can be calculated from those of half-cells. The forward reaction of Eq. 2 is spontaneous, which can be split into



Half-cell reaction	$E^{\circ}$ (V)
$\text{Zn}^{2+} + 2\text{e}^{-} \rightleftharpoons \text{Zn}$	-0.762
$\text{Fe}^{2+} + 2\text{e}^{-} \rightleftharpoons \text{Fe}$	-0.44
$2\text{H}^{+} + 2\text{e}^{-} \rightleftharpoons \text{H}_2$	0.0
$\text{Cu}^{2+} + 2\text{e}^{-} \rightleftharpoons \text{Cu}$	+0.34
$\text{Hg}_2\text{Cl}_2 + 2\text{e}^{-} \rightleftharpoons 2\text{Hg} + 2\text{Cl}^{-}$	+0.242

Select standard  
electrode potentials ( $E^{\circ}$ )

If,  $[\text{Cu}^{2+}] = [\text{Zn}^{2+}] = 0.1\text{M}$ ,

$$E = + 1.102 + (0.059/2) \times \log 0.1/ 0.1 = +1.102\text{V}$$

If,  $[\text{Cu}^{2+}] = 1\text{M}$  &  $[\text{Zn}^{2+}] = 0.1\text{M}$ ,

$$E = + 1.102 + (0.059/2) \times \log 1/ 0.1 = +1.13\text{V}$$

Typical Galvanic or Daniell cells generate  $\sim 1.1\text{V}$  unless the ion concentrations of copper & zinc are significantly different by a few orders of magnitude.

# **Questions and Discussion?**



# Biosensors

**B. Tech.**

**Course No.: EEL 3050**

**L-T-P [C]: 3-0-2 [4]**

***Prof. AJAY AGARWAL***


**ELECTRICAL ENGINEERING**

**IIT JODHPUR**

*Lecture 24 dated 17<sup>th</sup> Oct. 2024*

## Amperometric Biosensors:

- Enzymes are popular bioreceptors in identifying & quantifying their counterpart, substrates.
- Enzymes are typically proteins, while substrates are generally small chemicals.
- A good example is the use of glucose oxidase (GOx; the names of enzymes usually end with -ase) to quantifying glucose.
- GOx converts glucose into gluconic acid while generating hydrogen peroxide as a byproduct.
- During this reaction, GOx itself is reduced.



- An appropriate electron mediator can be added, to regenerate GOx back into its oxidized (i.e., active) form

- This **redox cycle generates** electrical **current** (with no voltage applied; not practical) or change in electrical current (with voltage applied).
- Since the **quantification** is made through **measuring current**, it is referred to as **amperometric** electrochemical **biosensor** or simply amperometric biosensor.
- This type of biosensor is particularly popular in measuring blood glucose level from diabetic patients.

## Conductometric Biosensors:

- In conductometric biosensors, conductivity (inverse of resistivity) is measured to monitor the redox reactions of enzymatic oxidation.
- It is similar to amperometric biosensors; the only difference is the current measured in amperometry & conductivity in conductometry.

$$\text{resistivity} = RA / l$$

where,

R resistance ( $\Omega$ )

A cross-sectional area of a specimen ( $\text{cm}^2$ )

l length of a specimen (cm)


The typical unit of resistivity is  $\Omega \cdot \text{cm}$ .

- Conductivity is defined as:  $\text{conductivity} = l / RA$


The unit of conductivity is  $\Omega^{-1} \cdot \text{cm}^{-1}$ , or  $\text{S} \cdot \text{cm}^{-1}$  (the standard unit), S, is called Siemens.

## Conductivity cell:


- Both resistivity & conductivity are measured with a **conductivity cell**,
- It has **two metal plates** separated by a **certain distance**.
- Here, **two  $1\text{ cm}^2$**  ( $= A$ ) metal plates separated by  **$1\text{ cm}$**  ( $= l$ ).
- The **ratio  $l/A$**  is referred to as **cell constant**, & it is  **$1\text{ cm}^{-1}$**  ( $= 1\text{ cm} / 1\text{ cm}^2$ )
- This **conductivity cell** is used to **quantify** the extent of the **enzymatic redox reactions**.
- Such **conductivity changes**, however, are generally **too small to be detected**.



- Conductivity is inverse of resistivity
- the small resistance change can be measured by a Wheatstone bridge
- If a conductivity cell replace the unknown resistor ( $R_4$ ) from a Wheatstone bridge, it is called conductivity bridge,
- a basic platform for many conductometric biosensors
- The conductivity cell is further miniaturized to provide better reproducibility & minimize the required sample volume.
- Interdigitated microelectrode (IDE) is one good example



Wheatstone Bridge



Conductivity Brid

# **Questions and Discussion?**