

# Biosensors

**B. Tech.**

**Course No.: EEL 3050**

**L-T-P [C]: 3-0-2 [4]**

***Prof. AJAY AGARWAL***


**ELECTRICAL ENGINEERING**

**IIT JODHPUR**

*Lecture 25 dated 21<sup>st</sup> Oct. 2024*

# Ion-Selective Electrodes (ISEs; Potentiometric):

- The **half-cells** can be used to measure the **concentration** of a **specific ion**.
- When a **half-cell** is used for ion concentration measurement, it is called an **ion-selective electrode (ISE)**.
- The **electrolyte solution** in the **electrode** makes contact with the **surrounding liquid** through a **membrane** that **allows** only a **specific type of ion** to pass through.



ISEs are classified in three types of membranes:

- i. solid-state,
- ii. liquid, and
- iii. Glass membrane.

i. **Solid-state ISEs** contain a crystalline membrane that is cut from a single crystal.


- Ex., a fluoride ISE uses a solid crystal of  $\text{LaF}_3$ , which allows only fluoride ions ( $\text{F}^-$ ) to pass through its membrane.
- The solid-state ISEs can also measure also include:  $\text{Ag}^+$ ,  $\text{Cl}^-$ ,  $\text{Br}^-$ ,  $\text{SCN}^-$ , and  $\text{S}^{2-}$ .

ii. **Liquid membrane ISEs** contain a plastic membrane, and the liquid ion-exchange material is absorbed into it.

- Vallinomycin-absorbed polyvinyl chloride (PVC) is an example of an ISE to selectively detect potassium ions ( $\text{K}^+$ ).
- Liquid membrane ISEs can be used for  $\text{NO}_3^-$ ,  $\text{Cu}^{2+}$ ,  $\text{Cl}^-$ ,  $\text{BF}_4^-$ ,  $\text{ClO}_4^-$ , &  $\text{K}^+$ .

iii. **Glass membrane ISEs** contain a **membrane** made from **thin glass** that is **very specific to hydrogen ions** ( $H^+$ ).

- The usual **composition** of the **glass** employed for detecting  $H^+$  is:  
22 %  $Na_2O$ , 6 %  $CaO$ , & 72 %  $SiO_2$ .
- **Glass membrane ISEs**, or **simply glass electrodes**, **can be** used to detect **other types of ions**, but they are primary used to measure  $H^+$ , or in other words, **pH**.
- The **detection limit** of **ISEs** is **very low**, ranging between  **$10^{-8}$  to  $10^{-11}$  M** (**10 nM to 10 pM**) of target ions.
- ISEs are **suitable** for measuring **low concentrations** in **small sample volumes**; since they do not chemically influence samples.



A glass electrode

## pH Electrode (Potentiometric):

- pH is a measure of  $H^+$  concentration.

$$pH = -\log a_{H^+} \approx -\log [H^+] \quad \dots \quad (17)$$


$$E = E^\circ + 0.059/n \times \log [M^+] \quad \dots \quad (9)$$

Plugging Eq. 17 into Eq. 9 gives

$$E = E^\circ - 0.059 \times pH$$


- That is **voltage decreases** by 0.059V or **59mV**, per unit increase in pH
- A **glass electrode** is used to **measure pH** if an **appropriate reference electrode** is used.
- A **calomel electrode** is most **frequently used** as such reference.
- The **two electrodes** can be **dipped** into a **soln.**, & the **voltage difference** between them can be measured, to **evaluate** the **pH** of the soln.

- the two electrodes are combined into a single electrode, known as a **combined pH electrode**
- At a **neutral pH of 7**, the combined pH electrode generates **0 mV**.
- If the **pH** of the solution **increases** by a **unit**, the **voltage** can **drop** by **59 mV** or vice versa
- the pH electrode produces 59 mV/pH unit, which is measured and displayed in pH units by the meter.
- The **meter** is a **voltmeter** that **displays** measurements in **pH units** instead of **volts**.
- pH electrodes generate a **lot of noise**, & **59 mV/pH** output is **low**, an **op-amp** is **required** in order to **construct a pH meter**.



A combined pH electrode

## pH meters:



Parts of a pH meter



# **Questions and Discussion?**