

Biosensors

B. Tech.

Course No.: EEL 3050

L-T-P [C]: 3-0-2 [4]

Prof. AJAY AGARWAL


ELECTRICAL ENGINEERING

IIT JODHPUR

Lecture 26 dated 23rd Oct. 2024


Lab-on-a-Chip Biosensors:

- A **lab-on-a-chip** (LOC) is a **device** that integrates **several laboratory functions** onto a **small platform**, typically only **millimeters** or **centimeters** in size.
- It is mainly a **network of channels & wells** that are **etched** onto **silicon** or **polymer** substrates to build **miniature laboratories**
- LOCs normally handle
 - **very small fluid volumes**, i.e. in **μL** or **nL** &
 - Are of **very small depth & width** of **channels**, **<1 mm**, i.e., in **μm** scale.
- **Pressure** or **electrokinetic forces** **move** small volumes of **liquid** in a **finely controlled manner** through the channels.




²
<http://commons.wikimedia.org/wiki/File:AutoFISH.jpg>


Microfluidics Platforms



Micro-Viscometer



Fluid FET



Lab-on-a-Chip (LOC) for Particles Separation & Counter

Ajay Agarwal and Balyan Prerna, Indian Patent No. 2592/DEL/2013, May 06, 2016

Filtration of Blood cells

- Separation of WBCs and RBCs from blood.
- Most of the WBCs were retained in the WBC separation and collection region, however some cells tends to escape due to deformable nature of the cells.





Figure: Separation of blood cells



Microfluidics is interdisciplinary

- related with various Science & Engineering areas



- Microfluidics deals with the manipulation/ control of fluids in μm scale.
- LOCs are also referred to as microfluidic devices or microfluidic chips.
- On a **single integrated system**, LOC enables sample
 - **handling, mixing, dilution,**
 - **electrophoresis & chromatographic separation,**
 - **staining, & detection**
- Electrophoresis separates components by their **rate of migration in an electric field**. It is often used to separate **macromolecules like proteins, DNA, & RNA** based on their charge, size, or binding affinity
- Chromatography separates components by **how they distribute** between a **stationary phase & a mobile phase**. The **mixture** is dissolved in a solvent & carried through a **system** with a **fixed material**. The components separate **because they have different affinities** for the stationary phase & are retained for **different lengths of time**.

Biosensors

B. Tech.

Course No.: EEL 3050

L-T-P [C]: 3-0-2 [4]

Prof. AJAY AGARWAL

ELECTRICAL ENGINEERING

IIT JODHPUR


Lecture 27 dated 24th Oct. 2024

Lab-on-a-Chip Biosensors:

The **main advantages** of the **LOC** are:

- ease of use,
- speed of analysis,
- low sample & reagent consumption, and
- high reproducibility due to **standardization & automation**.

Downscaling of volumes



Micro to Nano Fluidics

- When fluids flow in **nano-scale**, the flow is **non-continuous**
- Nanofluids are dispersions of **nanomaterials** in base fluids
 - nanoparticles, nanofibers, nanotubes, nanowires, nanorods, nanosheet, or droplets.
- Nanofluids have some **unique features** that are **quite different** from dispersions of mm or μm sized particles.
- Compared to conventional cooling liquids such as water, kerosene, ethylene glycol and microfluids, nanofluids have been shown to exhibit **higher thermal conductivities**.
- **Nanofluids do not block flow channels** and induces only a **very small pressure drop** during flow - beneficial for **heat transfer applications**

Applications of nanofluids:

(a) heat transfer,


(b) defect sensors,

(c) anti infection therapy,

(d) energy harvesting system ,


(e) hyperthermia and

(f) cosmetics



Microfluidic Devices fabrication

- MEMS technologies are used in realization of microfluidics' **key building blocks** such as **flow channels**, **pumps** and **valves**
- Typical microfluidic devices applications:
 - Chemical analysis, drug delivery, biological sensing, environmental monitoring, etc.
- Fluidic device **design considers** following effects:
 - type of flow (laminar or turbulent),
 - effect of bubbles,
 - capillary forces,
 - fluidic resistance &
 - capacitance



LOC fabrication with
photolithography

Biosensors

B. Tech.

Course No.: EEL 3050

L-T-P [C]: 3-0-2 [4]

Prof. AJAY AGARWAL


ELECTRICAL ENGINEERING

IIT JODHPUR


Lecture 28 dated 28th Oct. 2024

Microfluidic Devices

- i. **Flow channels:** A wide variety of microfluidic channels have been fabricated using
 - i. bulk micromachining (wet & dry etching),
 - ii. surface micromachining &
 - iii. moulding techniques



Soft lithography: replica molding



Various channels for microfluidic applications


ii. **Flow sensors:** MEMS flow sensors can be fluid-dependent flow or fluid-independent.

Fluid-dependent flow sensors measure the flow rate by heating a fluid 'upstream' & then recording its temperature 'downstream'.

- The flow rate is proportional to the temperature difference & transit time of the two actions.

Fluid-independent flow sensors measure pressure or force exerted on an object by the fluid.


- Using a **piezoresistive sensing mechanism**, flow measurements are fairly linear; direction and magnitude can be sensed by this method.
- In addition, fluid independent flow sensors do not involve any form of heating & hence are more suited to biological fluid applications.



An example of a MEMS bulk drag-force flow sensor

iii. Valves

- Valves are classified as **active** or **passive**; depending on **whether or not** they have an **external power** or **control source**. One of the simplest MEMS valves is the **passive** check valve



**Basic concept of
passive silicon
check valve**


- **Active** MEMS valves can be actuated by many methods: **thermal**, **piezoelectric**, **electrostatic** & **shape memory alloy** means. **Thermal actuation** is commonly used.

iv. Pumps

- Pumps are an important part of microfluidic devices
- MEMS pumps are very sensitive to fine particles which often cause contamination & leakage of the device
- Membrane, rotary & ultrasonic pumps are common types of MEMS pumps

Rotary pumps

- LIGA is used to fabricate micromachined PMMA gears in MEMS microfluidic systems
- These gears can be driven using electroplated NiFe bars mounted on the gears
- Fluid is pumped by the action of the turning gears




Basic concept of magnetic rotary micropump


vi. Droplet generators

- The majority of droplet generators in commercial MEMS microfluidic devices are inkjet printer heads
- These rely on either thermal or piezoelectric actuation to eject ink droplets;


- Thermal actuation are common
- Piezoelectrically actuated valves offer the advantage of very high forces, but very small movement for very large voltages
- They find use in automotive fuel-injection valves & spray nozzles.




The piezoelectric effect



An LOC made out of PDMS replica molding



Micro-transfer molding (μTM)



Micro-injection molding (μIM) and hot embossing

**Questions and
Discussion?**